

## FalisoMED Handicap Crank "LM"



|                 |          |
|-----------------|----------|
| Article number: | 03330255 |
|-----------------|----------|

|                  |            |
|------------------|------------|
| Original number: | 102 LM02 B |
|------------------|------------|

|      |           |
|------|-----------|
| EAN: | 999999991 |
|------|-----------|

|      |                 |
|------|-----------------|
| EVP: | <b>278.45 €</b> |
|------|-----------------|

### Technical data

|             |           |
|-------------|-----------|
| Main colour | silver    |
| Material 1  | aluminium |
| Weight      | 0.233 kg  |

### Description

#### Customize crank radius

The Crank 101/LM 02 (left) and Crank 102/LM 02 (right) crank extensions from FALISOMED can be used to individually adjust the range of motion of cyclists with restricted mobility or shorter legs.

The support block, which is mounted on the standard crank arm, can be individually positioned. For pedal arms of different diameters, there are holding blocks with different inner diameters. The retaining blocks can be individually positioned on the Crank Spider. The Crank Spider offers three different positions for the pedals.

The manufacturer recommends careful and precise assembly according to the instructions.

