

FalisoMED Handicap Crank "Reha"



Article number:	03330305
Original number:	101 REHA
EAN:	999999991
EVP:	424.90 €

Technical data

Main colour	silver
Material 1	aluminium
Weight	0.359 kg

Description

Auxiliary crank with pendulum arm

With the rehabilitation crank attachment Reha from FALISOMED, patients who have undergone surgery, for example, can improve their limited knee mobility by training on an Indoor Trainers or bicycle. The Crank Spider's range of motion is largely infinitely adjustable. The slide rail is available for the right (102 Reha) and left side (101 Reha) of the Crank Spider.

This means that the crank arm attachment is also suitable for use in physiotherapy practices or rehabilitation clinics.

The slide rail made of high-strength light metal is connected to the Crank Spider with two screw connections and a stainless steel dowel pin. A sliding block with a ball-bearing mounted pendulum crank is placed on the rail. The composition allows the pedal position to be adjusted practically infinitely and adapted to the training progress.

The manufacturer emphasizes the need for careful and precise assembly according to the instructions.