

## FalisoMED Handicap Crank "Reha-P"



|                 |          |
|-----------------|----------|
| Article number: | 03330339 |
|-----------------|----------|

|                  |           |
|------------------|-----------|
| Original number: | 102REHA P |
|------------------|-----------|

|      |           |
|------|-----------|
| EAN: | 999999991 |
|------|-----------|

|      |                 |
|------|-----------------|
| EVP: | <b>357.80 €</b> |
|------|-----------------|

### Technical data

|             |           |
|-------------|-----------|
| Main colour | silver    |
| Material 1  | aluminium |
| Weight      | 0.360 kg  |

### Description

#### Crank Spider for rehabilitation

With the rehabilitation crank attachments 101 REHA-P (for the left side) and 102 REHA P (for the right side) from FALISOMED, rehabilitation clinics or physiotherapeutic practices can respond to the limited hip or knee mobility of their patients, for example after an operation, and train their mobility again.

The Reha P crank consists of a Prisma sliding rail, which is connected to the crank without movement using two screw connections and a dowel pin. The sliding block sitting on the rail holds the pedals and can be adjusted in eight stages. This allows the respective Crank Spider radius to be adjusted accordingly.

The manufacturer recommends careful and precise assembly according to the instructions.